



# VIBRAGEM<sup>®</sup> DISCS

# BIOFREQUENCIES FOR CHILDREN AND ADULTS



# WHAT ARE VIBRAGEM<sup>®</sup> DISCS?

The discs are small, half-sphere shapes made of crystallized silicon (QUARTZ) and other minerals with a rigid, solid texture. During the manufacturing process, specific physical and environmental conditions, such as light, sound, and rest, were carefully considered. These conditions give VIBRAGEM<sup>®</sup> its unique properties, making it completely neutral.







**BIO-ENERGY RESONATOR**

**REGULATES ENERGY**


**BALANCES THE SYSTEM**

**RESTORES NATURAL HARMONY**

Vibrantum®








# HOW TO SENSE THE DISC'S VIBRATION?

- We will take deep breaths, focusing our attention on the area where the disc is placed.
- You may notice a slight tingling or a gentle warming sensation (each person experiences it differently).
- What you're feeling is your own energy being returned to you—amplified and regulated in frequency and wavelength—by the disc (this is one of the properties of quartz).

When unbalanced energy reaches the disc, it is returned amplified and regulated in both wavelength and frequency. This is one of the properties of quartz, which is why it is used in precision instruments.

Energy reaches the disc in various unbalanced wavelengths.



The disc acts as an amplifier and regulator of the body's waves and energy.

**Physical**

**25%**

**Mental**

**25%**



**25%**

**Emotional**

**25%**

**Energetic**



# APPLICATIONS



**PAIN**



**BEAUTY/BALANCE**



**DETOXIFICATION**



**REGENERATION**



**HARMONY**




**REJUVENATION**

# PHYSICAL ACTION

It provides more blood, oxygen, and resonance energy, so we can see its energetic and harmonious action through an infrared camera in the treatment of pain, stress, or imbalance. Its action on the natural balance of endorphins modifies the perception of pain or any physiological imbalance, playing an important role in repair and energy balance in the face of stress.


It has also been used with great effectiveness in the treatment of all types of back pain, migraines, sports injuries, body stress, and more.




# KNEE SURGERY



Before the therapy




6 hours of treatment






# LOCAL TREATMENT TO IMPROVE CIRCULATION AND FLUID ELIMINATION

The difference is noticeable after 24 hours of applying the disc, with temperature regulation, reduction of inflammation in the area, and improved circulation.



Before the therapy



24 hours of treatment








# Scar

We will place the discs on both ends of the scar for 7 days, then move them a little closer, and continue doing so until the scar disappears.






Recent scar from 2 months ago




Before the therapy



1 hour of treatment



24 hours of treatment



5 days of treatment

We place the discs on both sides to provide more energy.




Photo taken at 8 PM after a workday, showing inflammation caused by poor circulation.




Photo taken at 8 PM, 24 hours after applying the discs following a workday, showing a noticeable improvement in circulation.





# MEASUREMENT BEFORE AND AFTER THE APPLICATION OF THE DISCS.





A close-up photograph of a baby lying on their back, smiling and looking towards the camera. The baby has dark hair and eyes. They are wearing a white t-shirt with some faint, illegible text. On each of their upper arms, there is a small, square, white adhesive patch with a blue circular center. The baby is resting on a soft, textured surface, possibly a blanket or rug, with a warm, orange-toned background.

## **IT ACTS ON THE MORE SUBTLE ENERGIES.**

It is suitable for sensitive individuals,  
especially children and delicate skin.

We can create energetic geometric  
shapes.



# DISC THERAPY

- Painful areas and pressure points
- Programming and awareness
- Acupuncture protocols and bio-active points

- They are applied with hypoallergenic adhesive to painful points or selected areas until the pain, stress, stagnation, or energy blockage improves.
- This achieves constant stimulation and rebalancing in an easy manner. Therefore, knowledge of acupuncture points is not required to apply it. It also avoids the aversion that needles cause in some people, both children and adults.
- **In an acute crisis**, we can apply 5 to 6 discs during the first 24 hours, on acupuncture points and reflex zones, or in the area of acute pain. These discs are programmed according to the client's needs.






- **Permanent stimulation:** Wear two discs and place them on the abdominal area for 5 to 6 days, depending on the stress being treated. Renewal depends on the patient's emotional level and skin type. For example, when treating joint pain, we usually leave it on for 2 to 6 days, with an average of 3 days. For anxiety or stress, it's best to continue using it every day (as the problem improves, we can leave it on for fewer days).
- Leave it on for 4 to 6 days (it will not lose potency or effectiveness after 6 days; the change is made for hygiene reasons). Then, replace it with two new ones to maintain its vibrational effect.
- If using on young children, keep it out of their reach so they can't remove it.
- It can be used at any time: while sleeping, working, or meditating.
- It can be combined with other formulas if necessary.







The discs have no side effects. Additionally, being made of quartz, they have memory like water, which allows us to program them with information—whether through INTENTION, Laaps Wand, Color, Reiki, Flower Essences, Homeopathy, etc.

In this way, they open up an unlimited range of applications depending on the therapist's specialty.

# Vibrantum®



We are committed to providing solutions that not only improve the quality of life for those who use them but also promote a new way of caring for health. We would love to explore how Vibrantum® can become a strategic partner in your wellness product offerings.

[vibrantum638@gmail.com](mailto:vibrantum638@gmail.com)